

# PB 636



A special feature of the PB 636 bed is the possibility of lowering the mattress platform down to a level of 29 cm. The low position allows for increased comfort and significantly reduces the possibility of injury in the situation of accidental falls out. The height adjustment range allows caregivers to work in optimal conditions. The possibility of using the divided side rails means that the user does not have to be restricted by the side rails in the bed.



In the picture: PB 636 bed sonoma oak, with divided rails on one side of the bed. Accessories: bedside lamp.

# PB 636



PB 636 bed with divided railings on one side

## DETAILS:

1. Lowest mattress support platform level of 29 cm
2. Decorative cover of the mattress support platform
3. Latest actuator technology of bed lifting
4. Handset with lockable functions
5. Double brake castors
6. Divided side rails (option)

## OPTIONS:

length and mattress support material choices • colours choice • divided side rails

## BASIC TECHNICAL DATA:

Safe working load:	175,0 kg
Maximum user weight:	140,0 kg
Mattress support platform height adjustment	from 29 to 72 cm
Back rest adjustment:	0 ÷ 70°
Lower leg rest adjustment:	0 ÷ 20°
External dimensions (length × width):	207 × 106 cm
Clearance under the bed:	approx. 16 cm
Total weight:	118,6 kg
Castors:	100 mm with a brake


## FUNCTIONS:



ADJUSTABLE LEG  
SECTIONS



ADJUSTABLE BACK  
SECTIONS



ADJUSTABLE HEIGHT  
MATTRESS SUPPORT  
PLATFORM